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Nesrin BEN SALEM

Molecular Biology Product Manager Africa

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


1. Sample moves into lysis chamber. Cells and pathogens are lysed by bead beating, releasing nucleic acids
stranded DNA binding dye monitors each reaction

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Viruses

Not Detected Adenovirus F 40/41
Not Detected Astrovirus
Not Detected Norovirus GI/GII
Not Detected Rotavirus A
Not Detected Sapovirus

Run Details

Pouch: GI Panel v2.1
Run Status: Completed
Serial No.: 00788640
Lot No.: 133813

Protocol: Stool FA v2.3
Operator: John Madison (jrm)
Instrument: ITI FA "FA1315"

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Up to 22
Targets

• 14 bacteria
• 7 viruses

As target number varies between
any RP1.7, FilmArray RP2+, and
FilmArray RP EZ.

Blood Culture Identification Panel

27
Targets

- 19 bacteria
- 5 yeast
- 3 antibiotic-resistance genes

Gastrointestinal Panel

22
Targets

- 13 bacteria
- 5 viruses
- 4 parasites

Meningitis/Encephalitis Panel

14
Targets

- 6 bacteria
- 7 viruses
- 1 yeast

Pneumonia Panel

34
Targets

- 15 bacteria semi-quantitative
- 3 atypical bacteria
- 8 viruses + MERS-CoV
- 7 antibiotic-resistance genes

CE IVD and FDA-cleared

BIOFIRE FILMARRAY® RESPIRATORY 2 PLUS (RP2 PLUS) PANEL

BIOMÉ

FilmArray® Respiratory Panel 2 *plus* (RP2*plus*)

Fast

Accur

22
PATHOGENS
45^{min}

Viruses

Adenovirus
Coronavirus HKU1
Coronavirus NL63
Coronavirus 229E
Coronavirus OC43
Human Metapneumovirus
Human Rhinovirus/Enterovirus
Influenza A
Influenza A/H1
Influenza A/H3

Influenza A/H1-2009
Influenza B
Middle East Respiratory Syndrome
Coronavirus (MERS CoV)
Parainfluenza Virus 1
Parainfluenza Virus 2
Parainfluenza Virus 3
Parainfluenza Virus 4
Respiratory Syncytial Virus

Bacteria

Bordetella pertussis
Bordetella parapertussis
Chlamydia pneumoniae
Mycoplasma pneumoniae

BIOFIRE FILMARRAY® PNEUMONIA *PLUS* (PN *PLUS*) PANEL

BIOMÉ

BioFire® FilmArray® Pneumonia Panel

BACTERIA

(Semi-Quantitative)

Acinetobacter calcoaceticus-baumannii complex
Enterobacter cloacae complex
Escherichia coli
Haemophilus influenzae
Klebsiella aerogenes
Klebsiella oxytoca
Klebsiella pneumoniae group
Moraxella catarrhalis
Proteus spp.
Pseudomonas aeruginosa

ATYPICAL BACTERIA

(Qualitative)

Chlamydia pneumoniae
Legionella pneumophila
Mycoplasma pneumoniae

VIRUSES

Adenovirus
Coronavirus
Human Metapneumovirus
Human Rhinovirus/Enterovirus
Influenza A
Influenza B

ANTIMICROBIAL RESISTANCE

Methicillin Resistance
mecA/C and

Carbapenem

IMP
KPC
NDM
Oxa-48-like
VIM

ESBL

CTV M

Staphylococcus aureus
Streptococcus agalactiae
Streptococcus pneumoniae
Streptococcus pyogenes

Human Coronavirus
Syndrome Coronavirus
(MERS CoV)
Parainfluenza Virus
Respiratory Syncytial

Sample Type:

Sputum (including ETA) and BAL (including mini-BAL)

Accuracy

1 hr

- Enterococcus*
- Listeria monocytogenes*
- Staphylococcus*
 - Staphylococcus aureus*
- Streptococcus*
 - Streptococcus agalactiae*
 - Streptococcus pneumoniae*
 - Streptococcus pyogenes*

Acinetobacter baumannii
Haemophilus influenzae
Neisseria meningitidis
Pseudomonas aeruginosa
Enterobacteriaceae
Enterobacter cloacae complex
Escherichia coli
Klebsiella oxytoca
Klebsiella pneumoniae

Can
Can
Can
Can
Can
An
Re
med

BIOFIRE FILMARRAY GASTROINTESTINAL (GI) PANEL

BIOMÉ

FilmArray® Gastrointestinal (GI) Panel

Fast

Accurate

22

PATHOGENS

1 hr

Bacteria

Campylobacter
(*jejuni*, *coli*, and *upsaliensis*)
Clostridium difficile (toxin A/B)
Plesiomonas shigelloides
Salmonella
Yersinia enterocolitica
Vibrio
(*parahaemolyticus*, *vulnificus*
and *cholerae*)
Vibrio cholerae

Diarrheagenic *E. coli*/Shigella

Enteraggregative *E. coli* (EAEC)
Enteropathogenic *E. coli* (EPEC)
Enterotoxigenic *E. coli* (ETEC) *lt/st*
Shiga-like toxin-producing
E. coli (STEC) *stx1/stx2*
E. coli O157
Shigella/Enteroinvasive *E. coli* (EIEC)

Vi

Ad
As
No
Ro
Sa

Pa

Cr
Cy
En

BIO

BIOFIRE FILMARRAY® MENINGITIS / ENCEPHALITIS (ME) PANEL

FilmArray® Meningitis/Encephalitis (ME) Panel

Fast

Accurate

14
PATHOGENS
1 hr

Bacteria

Escherichia coli K1
Haemophilus influenzae
Listeria monocytogenes
Neisseria meningitidis
Streptococcus agalactiae
Streptococcus pneumoniae

Viruses

Cytomegalovirus (CMV)
Enterovirus
Herpes simplex virus 1 (HSV-1)
Herpes simplex virus 2 (HSV-2)
Human herpesvirus 6 (HHV-6)
Human parechovirus
Morbilli zoster virus (M/ZV)

Yeast

Cryptococcus neoformans

CE IVD

CE
IVD

⁷ Data on file at bioFire Diagnostics

CONCLUSION

Syndromic
approach

Why?

CONCLUSION

Syndromic
approach

Benefits

PATIENT

Provide better care faster
Avoid unnecessary use of
antimicrobials
Reduce length of hospitalization

CLINICIAN

Provide additional information
required to make informed
decisions on patient
management, often in critical
situations

COMMUNITY

Identifies infectious outbreaks
Helps reduce the inappropriate use
of antimicrobials and potentially
antibiotic resistance

PIONEERING DIAGNOSTICS